




WEIGHT GURUS



designer series
health metrics **scale**



Get Started: greatergoods.com/start



CUSTOMER SUPPORT

Email: info@greatergoods.com

Phone: (866) 991-8494

Website: greatergoods.com

Table of Contents

Important Safety Notes **2**

Introduction **3**

Scale Description **4**

Physical Features

Measuring Units

Setting the Measuring Unit

Display

Things To Know Before Using Your Scale **5**

Measuring

Auto-On

Auto-Detection

Press-Awake

User Setup **6**

SYNC With Your Scale **7**

Reading Your Results **8**

For Best Results

Typical Results **9**

Body weight

Body fat

Water weight

Muscle mass

Bone density

Troubleshooting **10**

Cleaning, Maintenance, and Disposal **11**

Technical Specifications **12**

Manufacturer's Warranty **13**

Important Safety Notes

Warnings

- Never use, or allow others to use this unit in combination with the following medical electronic devices:
 - Medical electronic implants such as pacemakers
 - Electronic life support systems such as an artificial heart/lung
 - Portable electronic medical devices such as an electrocardiograph
- This scale passes a harmless and unnoticeable electrical current through your body when taking a measurement. This electrical current is not felt while using the scale. This unit may cause the above mentioned medical devices to malfunction.
- This product is not intended for use by pregnant women.
- This product is not intended for use by infants, toddlers and children under 10 years of age.
- Do not step on the edge of the scale while getting on or off, otherwise it may tip.
- Do not jump on scale.
- Protect scale from hard knocks, temperature fluctuations and heat sources that are too close (e.g. Stoves, heating units)
- Do not drop scale or drop any objects on it as this may damage the sensors.
- Do not step or stand on the scale when your body and/or feet are wet. For example, after taking a shower.
- Do not step on the scale when the surface is wet. There is a danger of slipping.
- Please be aware that the measurements obtained from this body analysis scale represent only an approximation.
- Do not step on this scale if you weigh more than 400 lbs.
- Place the scale on level flooring, otherwise it may tip.
- Clean the scale with a slightly damp cloth, DO NOT immerse the scale in water or use chemical/abrasive cleaning agents.

Introduction

Welcome to your new Weight Gurus Scale.

Being healthy is worth the effort! That's why we've designed the Weight Gurus Scale to be more than just another scale. It's a beautifully designed, fully-featured body composition scale that, when desired, connects to the power and convenience of your smart phone. This lets you track your progress over time, helping you stay motivated and reach your goals.


To automatically sync your results to your smart phone, download the free Weight GURUS App for iOS (Apple) or Android (see more information on page 7).



For more product information, please visit GreaterGoods.com

Scale Description

Physical Features



Measuring Units


lb	pound
kg	kilogram

Setting the measuring unit

By pressing the button on the back of the scale, you can switch between lb. (pound) and kg (kilogram).

Display

- SYNC
- Body Fat
- Muscle Mass
- Water Weight
- Bone Density



Things To Know Before Using Your Scale

Measuring

Measuring your weight and calculating body composition is easy. Simply step on your scale with bare feet, and evenly distribute your weight. Make sure that your feet are making good contact with all four metal conductors. Body composition can only be displayed if you have properly set up a user profile on your scale (page 6).

Auto-On

Auto-on is the easiest way to use your scale. Simply step on the weighing surface with both feet and evenly distribute your weight. Once you complete user setup (page 6), Auto-Detection will match you to your profile automatically. If no match is made to a user profile, then the scale will only display your weight information. Body composition calculations require specific variables such as your height, gender and activity level, so a user profile must be setup in order to display those measurements.

Auto-Detection

Once a profile has been setup, your scale will automatically match you to your correct profile by weight. The scale determines who the user is by a range of 10 lbs. and in some cases a conflict may arise if two or more users are within ten lbs. of each other (see "Troubleshooting" for resolving an auto-detection conflict).

Press-Awake

Another easy way to turn on your scale is to briefly step-on, then off the weighing area of your scale. "Press-Awake" allows you to enter setup mode to create or edit a new user (see User Setup, page 6). If needed, it also allows you to select your user profile manually.

User Setup

IMPORTANT: Place your scale on a hard flat surface, and be ready with bare feet to record your first weight entry on the last step.

This scale can be customized for up to 8 users.

Scale Setup

1. Press-Awake your scale on. Simply step on then off to begin.
2. Press and hold SET.
3. Use ▲ and ▼ to make selections, then tap SET to select and continue for each option below.

User Number	Choose a user profile number, tap SET
Gender	Choose Male or Female, tap SET
Athlete Mode	Choose Athlete Y or N, tap SET (See Athlete Mode description below)
Height	Select your height, tap SET
Age	Select your age, tap SET
IMPORTANT: YOU MUST NOW WEIGH TO RECORD YOUR INITIAL WEIGHT ENTRY	
Weight	Step on the scale with bare feet, evenly distributing your weight. Make sure that your feet are making good contact with both metal conductors.


*At any time in the setup process, if you have to leave or the screen shuts off due to inactivity, just start over from number one above.

Athlete Mode

An athlete is defined as a person who is involved in intense physical activity approximately 12 hours per week and who has a resting heart rate of approximately 60 beats or less per minute.

Sync With Your Scale

Syncing your results to your smart phone is easy.



1. Get the App
2. Step on the scale (You must use an active user profile to SYNC Body Composition, see page 6 for setup instructions)
3. After the scale completes measuring, press the SYNC button to display your measurements in barcode format.
4. Use the Weight GURUS App and your smart phone's camera to scan the barcode. Make sure to align the brackets on your phone's screen to the scale's barcode screen.

TIPS:

ZOOM: Use the zoom buttons on the left to get closer or further away to properly align the app screen's brackets with the scale's screen.

GLARE: Lean from side to side when scanning to avoid direct glare from overhead lights. If a bright reflection is blocking your view of the barcode, the scanner may not detect your results.

*With the Weight GURUS App, Weight, Body Fat % & Muscle Mass % are automatically transferred to the App via bar-code display.

Learn more at GreaterGoods.com

Reading Your Results

Once you have properly weighed yourself, the scale will cycle through the measurements it has calculated for you.

Weight	Body Fat	Muscle Mass	Water Weight	Bone Mass
126.2	20.7	38.2	56.9	4.0

For best results

Remember to place your scale on a hard level surface. Step on the scale with bare feet, evenly distributing your weight. Make sure that your bare feet are making good contact with all four metal conductors.

Using this scale on carpet will not give accurate results.

Wearing socks during weighing will result in incorrect body composition measurements.

It's important to follow a routine. We recommend weighing daily around the same time each day. Alternatively, you can weigh weekly. When weighing weekly we recommend weighing on the same day each week.

Weighing under the same conditions will help ensure accurate results. Weight can vary due to different clothing selections, weighing at different times of day, weighing before or after meals, etc. Following a routine will help achieve accurate results.

Typical Results

Keep in mind that you know your body best. The following ranges of Body Fat %, Water Weight %, Muscle Mass % and Bone Mass % are offered only as guides. The best way to understand your body composition is to watch for changes. For example, watch to see how your body fat is going up or down over time.

Body Fat %

Age	Female				Male			
	Low	Normal	Above Avg.	High	Low	Normal	Above Avg.	High
20-29	<19	20-28	29-31	>31	<13	14-20	21-23	>23
30-39	<20	21-29	30-32	>32	<14	15-21	22-24	>24
40-49	<21	22-30	31-33	>33	<16	17-23	24-26	>26
50-59	<22	23-31	32-34	>34	<17	18-24	25-27	>27
60+	<23	24-32	33-35	>35	<18	19-25	26-28	>28

Water Weight %

Age	Body Fat %	Water %
Male	4 - 14%	63 - 70%
	15 - 21%	57 - 63%
	22 - 24%	55 - 57%
	25 and over	37 - 55%
Female	4 - 20%	58 - 70%
	21 - 29%	52 - 58%
	30 - 32%	49 - 52%
	33 and over	37 - 49%

Muscle Mass %

Female	Approx. 30%
Male	Approx. 40%

Bone Mass %

	Female			Male		
Weight	<100	100-135	>135	<135	135-165	>165
Bone	3.9	4.1	4	4	4.4	4.2

Troubleshooting

Low battery

“Lo” is displayed when the batteries need to be replaced. Open the battery door on the back of the scale, and remove and replace the batteries with four AAA batteries.

Overload

Overload occurs when the scale has too much weight (Capacity: 180kg / 397lb). Please remove the weight on the scale to protect the scale sensors.

Measuring error

If the scale has trouble measuring your weight, it will show “Err” and turn off. This is most likely due to not settling on the scale fast enough. Make sure that the scale is placed on a flat hard surface (not carpet), and you distribute your weight evenly while holding still.

Failed to auto-detect your user profile while weighing

The weight data is shown only and the scale will turn off after a few seconds. Please manually select your number and measure once again (Press-Awake, select user number, Tap SEL).

Resolving an auto-detection conflict

An auto-detection conflict will occur if two or more users are within 10lbs. of each other. When an auto-detection conflict occurs, your scale will give you a chance to select your profile number manually. Use the arrows to select your profile, tap SEL, and step on the scale. The scale will now proceed measuring for the profile you have selected.

How do I clear a user profile?

First, Press-Awake the scale on, then press and hold SEL. Now, select the user to delete and press SEL. Once selected, press the up arrow to change Clr to “Y.” Now press SEL to clear the user. The user is now removed.

Cleaning, Maintenance and Disposal

Cleaning

Your scale is not waterproof. Never immerse the scale in water or place it in a dishwasher. To clean the scale platform, use a soft, slightly damp cloth or sponge. For stains or other residue, use a mild dish washing soap. Never use harsh cleaners or treated cloths.

Changing the batteries

Your scale uses 4 AAA batteries. To change the batteries, open the battery door on the bottom of the scale, remove the old batteries and replace them with the + end oriented as shown in the battery compartment.

Disposal

Batteries should be disposed of separately from household waste. Always dispose of batteries as per your local regulations. Dispose of this device in accordance with your local regulations.

Technical Specifications

Weight unit: lb / kg

Capacity: 397lb / 180kg

Division: 0.1lb / 0.1kg

Minimal weight: 5kg

Function keys: SYNC, ▲, SEL, ▼

Display: Positive LCD

Batteries: 4 AAA

User No.: 8

Age range: 10 - 85

Stature range: 3ft 3inches - 7ft 3inches / 100-220cm

Measurements calculate: Body weight, body fat, water weight, muscle mass and bone mass

Working temperature: 32°F-104°F / 0°C-40°C

working humidity: ≤90% RH

Storage temperature: -4°F - 140°F / -20°C - 60°C

Manufacturer's Warranty

Your scale is warranted by the manufacturer against defects in materials and workmanship for five (5) years from the original purchaser from the date of purchase. Proof of purchase is required.

The warranty is void if the product has been subjected to mechanical damage or mistreatment, such as immersion. This warranty is in lieu of all other warranties, and limits the liability of the manufacturer. This warranty gives you certain legal rights and you may have other rights depending on which state the product was purchased.

If your scale is defective, please contact **Greater Goods, LLC**.

Phone: (866) 991-8494


Email: info@greatergoods.com

Web: GreaterGoods.com

Greater Goods, LLC
125 N. Main St.
Suite 202
St. Charles, MO 63301



QUICK START GUIDE




1. SETTING UP YOUR SCALE

- A** Remove the plastic pull tab from back of scale.
- B** Tap-Awake the scale on. (Briefly step on, then off the scale.)
- C** Press and hold SEL.
- D** Use ▲ & ▼ to make selections, then tap SEL to select.
Enter: USER NUMBER*, GENDER, ATHLETE MODE*, HEIGHT, AGE.

* **USER NUMBER** - Your personal user profile number (Select user 1-8)
ATHLETE MODE - An athlete is defined as a person who is involved in intense physical activity approximately 12 hours per week and who has a resting heart rate of approximately 60 beats or less per minute.

2. WEIGH

Step on the scale with bare feet, evenly distributing your weight. Make sure that your feet are making good contact with all four metal conductors.



3. SYNC

- A** Press SYNC on scale as it's showing measurement results.
- B** Open the Weight GURUS App's Camera/Sync function.
- C** Align the brackets on your phone's screen with the scale's LCD.
When the display shows the barcode, results will be automatically transferred to your smartphone. Press ENTER to store your results.

For more info visit: www.greatergoods.com